


MODEL 9302 & MODEL 9304



Model M9302 & Model M9304 Temperature Indicators

Two or four input channels with Ethernet connectivity Scientific Instruments' Model 9302 and 9304 are the most flexible and accurate temperature monitors currently available. Virtually any cryogenic temperature sensor from any manufacturer can be selected by a single setting of the front panel. Additional custom or specially calibrated sensors require only a simple setup procedure. In addition to its high accuracy/performance and low noise design, unique features include: Constant Voltage AC sensor excitation, Internal Data Logging, Ethernet Connectivity, a large easy to read display and extensive utility software.

Ethernet connectivity adds a new dimension of utility to these monitors. In both industrial and laboratory applications, Ethernet is more reliable and easier to use compared to other communication standards. Furthermore, it is essential to remote, distributed sensor or Local Area Network based systems.

- Two models available. The Model 9304 has four input channels and Model 9302 has two. All channels are identical in function.
- Multipurpose input channels support Silicon Diode, Ga, Platinum RTD and Cryogenic NTC resistive temperature sensors.
- Constant-Voltage, AC excitation of resistive sensors increases temperature range and improves sensitivity.
- Analog voltage output plus two programmable dry-contact relays.
- Continuous data logging into internal Non-Volatile Memory.
- High speed Ethernet interface. Electrically isolated.
- Built-in web server. Temperature monitoring and instrument configuration can be performed using any web browser.
- TCP/IP User Data Socket for complete remote operation using a simple IEEE-488 like command language.
- SMTP interface sends e-mail on a selected alarm condition.
- LabView drivers available for both the Ethernet (TCP/IP) and serial port.

Sensor Inputs: The Model 9302 & 9304 can be easily configured to support virtually any cryogenic temperature sensor. Configuration is performed from the instrument's front panel or a remote interface. There are no jumpers, trim pots or switches.

Silicon Diode sensors from Scientific Instruments or any other manufacturer are directly supported over their full 1.4 to 500K range using built-in calibration curves and sensor data. Plus, non-volatile Flash memory is available for several custom or calibrated sensors.

Platinum RTD sensors can use built-in DIN 43760 (IEC 750) standard setups for 100Ω OR 1000Ω devices. The Model 9302 & 9304 uses the DIN standard for temperatures from 70K to 1020K and extends it down to 30K for cryogenic use. Operation down to about 14K is available using user supplied curves.

A unique feature of these monitors is the use of a ratiometric resistance bridge technique to measure Platinum RTD sensors. This significantly reduces low frequency noise and drift to provide solid measurements.

These monitors provide robust support for the Negative Temperature Coefficient (NTC) sensors commonly used by cryogenic applications. They include *Ruthenium-oxide*, *Cernox™*, *Carbon-Glass™*, *Germanium* and several others. Since they have a negative temperature coefficient, the constant-voltage measurement method will reduce, rather than increase, power dissipation in the sensor as temperature decreases. By maintaining the lowest possible power level, sensor self-heating is minimized and useful temperature range is greatly increased.

An additional advantage to constant-voltage biasing is that NTC resistors lose sensitivity in the upper part of their range. By auto-ranging excitation current to maintain a constant voltage, sensitivity and noise immunity in that range is greatly improved.

Sensor excitation used in conjunction with the constant-voltage feature is a 2.5 Hz bipolar square wave. This effectively cancels thermal EMF induced offset errors that sometimes occur in cryogenic systems. The maximum and minimum sensor resistance that can be read is a function of the selected voltage bias.

Resistance Range Table		
Voltage Bias	Min Resistance	Max. Resistance
10.0mV	10Ω	1.0MΩ
3.33mV	3.3Ω	430KΩ
1.0mV	1Ω	100KΩ

Accuracy: Measurement accuracy is obtained by using a 24-bit analog to digital conversion. Accuracy is further enhanced by extensive use of Digital Signal Processing (DSP) techniques.

The Model 9302 & 9304 include built-in curves that support most industry standard temperature sensors. Additionally, four user calibration curves are available for custom or calibrated sensors. Each curve may have up to 200 entries. The accuracy of any sensor can be greatly improved by the use of this feature. This will fit a Diode, Platinum RTD or NTC resistor sensor's calibration curve at up to three user specified temperature points. It is a built-in, easy to use method for obtaining higher accuracy temperature measurements without expensive sensor calibrations.

Lowest Noise: These temperature indicators were designed for use in extremely low noise environments that cryogenic systems often require. Linear-mode power supplies are used throughout and sensor excitation current sources are not multiplexed. The enclosure is all Aluminum with wide conductive overlaps on all mating metal surfaces so that radiated RFI noise is virtually eliminated. An effective shielding and grounding scheme further allows the user to minimize both conducted and radiated noise.

Easy to use: The Model 9302 & 9304 front panels consist of a large, bright Vacuum Fluorescent display and a 5-key keypad. Most features and functions can be accessed via this simple and intuitive menu driven interface. Two temperature readings can be displayed in a large, easy to read 10mm two-line font. Underlying menus switch to a 5mm high four-line font for more content. Additionally, the Model 9304 can display all four inputs plus alarm and relay status in this smaller font. Temperature displays are autoranged to show the most number of significant digits. Built-in digital filters can be used to smooth temperature data. Displays are in units of K, °C, °F, Volts, or Ohms. The status of built-in alarms and relays are indicated by LEDs located to the right of the display.

Outputs: The Model 9302 & 9304 each have two dry-contact relay outputs. Either may be independently programmed to assert or clear based on a high or low temperature condition. Normally-open contacts are available on the rear panel. Also available is a single analog output channel. This is a zero to 4 Volt output that is proportional to any selected input.

Data logging is performed by continuously recording temperature and status to an internal 20K byte circular memory buffer. Data is time stamped so that the actual time of an event can be determined. Non-volatile memory is used so that data will survive a power failure.

Ethernet: The Model 9302 & 9304 connects directly to any 10-Based T Ethernet interface to make measurements easily and economically. Simple connection to any existing Local Area Network allows stable, precise, cost-effective measurements in laboratory or industrial environments as well as in remote, distributed data acquisition systems. The Model 9302 and 9304 can even be connected directly to the Internet with a user-supplied IP address.

Using Ethernet HTTP protocol, the monitor's embedded web server allows the instrument to be viewed and configured from any web browser.

Input channels can be configured using text entry and drop-down box selections by going to the monitor's "inputs" web page.

Using SMTP protocol, the monitor will send e-mail based on selected alarm conditions. E-mail is completely configured by using the web page interface.

The TCP/IP data port server brings fast Ethernet connectivity to all common data acquisition software programs including LabView™.

The TCP/IP protocol is used to implement a text based command language like those commonly used with IEEE-488 or RS-232 interfaces. This is the primary way that user software interfaces to the monitor.

The remote command language is SCPI compliant according to the IEEE specification.

With Ethernet connectivity, the user has complete control of the monitor by using any web-enabled device from desktop PC to a wireless Pocket PC™. It is platform and operating system independent, working equally well with Windows, Linux or Macintosh based computers. There are no expensive cards or cables and, no confusing configuration requirements.


Software: Utility software is provided that connects any Windows based personal computer to either monitor. This software provides a graphical control panel that greatly simplifies instrument setup and configuration.

Sensor calibration curves may be downloaded to the monitor, viewed and edited.

LabView: LabView™ drivers are provided for both the Ethernet TCP/IP and RS-232 interfaces.

Rear Panel Connections

- Input Connectors: Two DB-9 receptacles provide 4-wire measurement connection to two sensors each.
- LAN: Standard RJ-45 Ethernet connector with built-in connection and activity LEDs.
- RS-232: Null-Modem connector (DB-9, pins).
- Relays / Analog Output: 6-pin detachable terminal block 3.5mm.
- AC Power: RFI filtered Power Entry Module including AC power line switch and fuse drawer. Line voltage selection is performed by Internal jumpers.



Sensor Performance

Silicone Diode sensors use a fixed excitation current of 10 μ A/100 μ A and an input voltage range of 0 to 2.5V.

Both Negative and Positive Temperature Coefficient (NTC) / (PTC) resistor sensors are supported using a ratiometric bridge technique to cancel low frequency noise.

Other sensor types include: Platinum, Excitation currents are 1.0mA, 100 μ A and 10 μ ADC. Corresponding full-scale resistance ranges are: 312 Ω , 3.2k Ω and 31k Ω .

NTC sensors include: Ruthenium Oxide and Cernox™, Constant-voltage AC sensor excitation allows the use of these sensors over an extended temperature range. Excitation voltage selections are 10mV and 1.0mV.

Sensor Performance Data					
Sensor Type	Silicone Diode	100 Ω Platinum DIN43760	1000 Ω Platinum DIN43760	Ruthenium Oxide1	Cernox™ 1
Sensor Sensitivity	300K: 2.4mV/K 77K: 1.9mV/K 4.2K: 30mV/K	800K: 0.36 Ω /K 300K: 0.39 Ω /K 77K: 0.42 Ω /K 30K: 0.19 Ω /K	600K: 3.7 Ω /K 300K: 3.9 Ω /K 77K: 4.2 Ω /K 30K: 1.9 Ω /K	1.0K: 1260 Ω /K 4.2K: 80.3 Ω /K 20K: 3.96 Ω /K	1.4K: 240K Ω /K 4.2K: 2290 Ω /K 77K: 2.15 Ω /K 300K: 0.16 Ω /K
Measurements Accuracy	300K: 21 μ V 77K: 23 μ V 4.2K: 44 μ V	800K: 2.4m Ω 300K: 2.4m Ω 77K: 1.2m Ω 30K: 1.2m Ω	600K: 38m Ω 300K: 38m Ω 77K: 4.7m Ω 30K: 4.7m Ω	1.0K: 1.9 Ω 4.2K: 1.4 Ω 20K: 1.09 Ω	1.4K: 675 Ω 4.2K: 5.1 Ω 77K: 161m Ω 300K: 450m Ω
Temperature Measurement Accuracy	300K: 8.7m 77K: 12mK 4.2K: 1.6mK	800K: 6.7mK 300K: 6.2mK 77K: 2.8mK 30K: 9.8mK	600K: 6.2mK 300K: 6.2mK 77K: 2.8mK 30K: 9.8mK	1.0K: 1.9mK 4.2K: 17mK 20K: 275mK	1.4K: 2.2mK 4.2K: 2.2mK 77K: 75mK 300K: 295mK
Measurement Resolution	300K: 7.4 μ V 77K: 7.4 μ V 4.2K: 15 μ K	800K: 1.8m Ω 300K: 1.8m Ω 77K: 460m Ω 30K: 460 μ Ω	600K: 15m Ω 300K: 15m Ω 77K: 1.8m Ω 30K: 1.8m Ω	2.0K: 11m 4.2K: 11m Ω 20K: 11m Ω	4.2K: 46m Ω 77K: 1.8m Ω 300K: 0.5m Ω
Temperature Resolution	300K: 3.0mK 77K: 3.8mK 4.2K: 500 μ K	800K: 1.8m Ω 300K: 4.7mK 77K: 1.1mK 30K: 2.4mK	600K: 4mK 300K: 4mK 77K: 0.5mK 30K: 1.0mK	2.0K: 30 μ K 4.2K: 0.13mK 20K: 2.9mK	4.2K: 50 μ K 77K: 0.85mK 300K: 3.5mK
Power Dissipation	4.2K: 17 μ W 77K: 12 μ W	30K: 3.7 Ω W 77K: 20 μ W	30K: 370nW 77K: 20 μ W	1.0K: 42nW 4.2K: 73nW	1.4K: 1.1nW 4.2K: 20nW
Magneto-resistance	Very Large	Moderate	Moderate	<2% for H<2T	<1% for H<2T

1 10mV Constant-Voltage excitation.

Specifications

USER INTERFACE

Display Type: Graphics VFD, 10mm character height.

Number of Inputs Displayed: Two or Four.

Keypad: Sealed Silicon Rubber.

Temperature Display: Six significant digits, autoranged.

Display Update Rate: 0.5 Seconds.

Display Units: K, C, F or native sensor units.

Display Resolution: User selectable to seven significant digits.

INPUT CHANNELS

There are two input channels on the Model 9302 and four on the Model 9304, each may be independently configured for any of the supported sensor types.

Sensor Connection: 4-wire differential. DB-9 receptacle.

Sensor Types: See Supported Sensor Table.

Sensor Selection: Front Panel of remote interface.

Sensor Resolution: Sensor Dependent. See Sensor Performance Data Table.

Sensor Excitation: Constant current: 1mA, 100 μ A or 10 μ A.

Constant Voltage: 10mV, 3.3mV, and 1.0mV RMS with excitation currents from 1.0mA to 10nA in steps of 5% of power.

Resistance Measurement type: Ratiometric bridge.

Resistance Range: Constant-voltage resistance measurement range 10 Ω to 350K Ω .

AC Excitation Frequency: Resistor sensors in constant-voltage mode: 2.5Hz bipolar square wave.

Sample Rate: 10Hz per channel.

Measurement Resolution: Sensor dependent.

See Sensor Performance Data Sheet.

Digital Resolution: 24 bits.

Digital Accuracy: 0.0015% of full scale.

Measurement Drift: <15ppm/°C.

Measurement Filter: 0.5, 1, 2, 4, 8, 16, 32, and 64 Seconds.

Calibration Curves: Built-in curves for industry standard sensors plus four user curves with up to 200 entries each. Interpolation is performed using as Cubic Spline.

DATA LOGGING

Data logging is performed to an internal, 20K-byte circular buffer and is time-stamped with a real-time clock.

Buffer memory is non-volatile and will retain valid data without AC power.

ANALOG OUTPUT

The analog output is a scaled voltage output that is proportional to any selected input temperature.

Output Range: Zero to 4 Volts.

Output Impedance: 500 Ohms.

Digital Resolution: 0.0015% of full-scale range.

Connection: Detachable terminal block.

RELAY OUTPUTS

Each relay output may be programmed to assert upon detection of a high or low temperature on any selected input channel.

Number: Two. Dry, Normally Open contacts.

Contact Rating: 30VDC at 1A. Connector.

Connection: 6-pin detachable terminal block.

REMOTE INTERFACES

Ethernet: 10-Base T. Electrically isolated

TCP/IP user data socket provides remote control and data interface to common data acquisition software by using an ASCII command language
HTTP provides built-in web server.

SMTP sends e-mail based on user selected alarm conditions.

RS-232: Serial port as an RS-232 standard null modem.

Data Rates are 9600, 38,400 and 57,600 Baud.

Language: Remote interface is IEEE SCPI compliant on both the TCP/IP and RS-232 Interface.

LabView® Drivers are available for the Ethernet TCP Data Socket and the RS-232 interface.

GENERAL

Ambient Temperature: 25°C \pm 5°C for specified accuracy

Mechanical: 8.5"W x 3.5"H x 12"D. One half-width 2U rack.

Instrument bail standard, rack mount kit optional.

Weight: 5 Lbs.

AC Power Requirement: 110 of 220VAC, +5% to -10%, 50 to 60Hz, 30VA.

AC Power Switch: Rear Panel

Conformity: European CE certified

Options & Accessories

Included Accessories

- User's Manual
- Utility Software CD
- Relay Connector
- Dual sensor input connector/cable assembly
- Detachable 120VAC Line Cord
- Certificate of Calibration

Optional Accessories

- Panel Mount hardware
- Null Modem Cable
- Shielded sensor connector kit.

Ordering Information

Part Number	Description
9302-110	Two-channel monitor set for 90 to 120VAC
9302-240	Two-channel monitor set for 220 to 240VAC
9304-110	Four-channel monitor set for 90 to 120VAC
9304-240	Four-channel monitor set for 220 to 240VAC



4400 West Tiffany Drive, West Palm Beach, Florida 33407
 Telephone: (561) 881-8500 • Fax: (561) 881-8556
 Visit our website @ www.scientificinstruments.com
 e-mail: sales@scientificinstruments.com