

Sensor Type	Thermocouple
Model Number	Chromel-Gold
Serial Number	(Standard)
Current	10 μ A
Gain	1
Temperature	Data
1	-5.30145
2	-5.29203
3	-5.28126
4	-5.26934
5	-5.25644
6	-5.24271
7	-5.22827
8	-5.21326
9	-5.19778
10	-5.1819
11	-5.16571
12	-5.14927
13	-5.13265
14	-5.11588
15	-5.09901
16	-5.08207
17	-5.06509
18	-5.0481
19	-5.03111
20	-5.01413
21	-4.99718
22	-4.98026
23	-4.96338
24	-4.94655
25	-4.92976
26	-4.91302
27	-4.89632
28	-4.87966
29	-4.86304
30	-4.84646
32	-4.81338
34	-4.7804
36	-4.74747
38	-4.71457

40	-4.68164
42	-4.64867
44	-4.61563
46	-4.58249
48	-4.54924
50	-4.51585
55	-4.43172
60	-4.34656
65	-4.26031
70	-4.17298
80	-3.99508
90	-3.81298
100	-3.62684
110	-3.43687
120	-3.24339
130	-3.04672
140	-2.84715
150	-2.64491
160	-2.44018
170	-2.23313
180	-2.02395
190	-1.81281
200	-1.59985
210	-1.38518
220	-1.16894
230	-0.95129
240	-0.73249
250	-0.51272
260	-0.29197
270	-0.07011
273.15	0
280	0.15264
290	0.3755
300	0.5992
400	2.79825
446	3.8729
500	5.1347
600	7.4697